

Mass spectrometry applying soft photo-ionisation for real time characterisation of transients from flash pyrolysis of biomass

Thorsten Streibel


University of Rostock, Institute of Chemistry, Analytical Chemistry

Workshop at 19th EU Biomass Conference
ICC Berlin, 08.06. 2011

Intention and purpose for fast MS analysis in biomass pyrolysis

- Which species are present in pyrolysis gas prior to condensation?
- Does the composition of pyrolysis gas influence product quality?
- Is it feasible to employ real time methods for characterisation and even control of biomass pyrolysis?
- Selected application: Technical biomass pyrolysis facility at Karlsruhe Institute of Technology (first step of bioliq[®] process)
- Monitoring of gaseous pyrolysis products of a large variety of biomass (wood, straw, rapeseed residue, corn cob, Miscanthus, palm frond etc.)


Experimental setup



Scheme of the on-line process analysis

Ionisation by absorption of two **UV-photons**

⊗ Resonance Enhanced Multiphoton Ionisation (REMPI)



- Photon-Energy 4 - 6 eV

- **Selective and sensitive ionisation of (poly)aromatic compounds**

Ionisation by absorption of one **VUV-photon**

⊕ Single photon ionisation (SPI)

2X ————— $M^+ + e^-$ Single Photon Ionisation:



- Photon energy 8 - 11 eV

- „General“ ionisation of organic compounds

VUV- und UV- photon sources

Laser generated UV-photons

- Nd:YAG Laser (266 nm)

- compact
- high energy density

Pulsed ionisation

Lamp generated VUV-photons

- Electron beam pumped rare gas excimer lamp (EBEL)


- brilliant light source
- high photon yield
- robust and easy to operate

Continuous ionisation

Flash pyrolysis of wheat straw




REMPI-MS




SPI-MS

● Ketons/Aldehydes






Flash pyrolysis of rapeseed residue



Current status and future work

- On-line PI-MS successfully adapted and applied for detection of gaseous biomass pyrolysis products
- Correlation of on-line data with respective composition of bio-oil and/or tar
- Analysis of these condensation products imaginable by
 - Pyrolysis combined with GCxPI-MS or GC-MS or PI-MS
 - Comprehensive two-dimensional GCxGC
 - High resolution mass spectrometry





Second column

